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INSTRUCTIONS FOR USE KNEE IMMOBILIZER "IMMOBIMAX" O-1250

CAUTION! Before using the device, please consult the instructions for use.

I. DEVICE DESCRIPTION


The "IMMOBIMAX" O-1250 knee immobilizer consists of a body made of soft, perforated material, four aluminum rails (two in the rear part of the brace and two located on the outside, far apart from each other) and straps. The aluminum rails are placed in pockets, so that they can be removed. Thanks to the aluminum rails, the "IMMOBIMAX" O-1250 brace can be easily adjusted to the shape of the leg, which ensures proper stabilization and comfort of use.

II. FEATURES


The "IMMOBIMAX" O-1250 knee immobilizer, thanks to its rigid and open construction, is a good solution for use after various types of procedures in the knee joint. "IMMOBIMAX" O-1250 enables very fast knee stiffening. The application is easy enough for anyone to handle it, in the least expected situation. Thanks to the stiffening rails, it can be used, among others, during postoperative rehabilitation, after injury to the patella or patellar ligament as well as during the treatment of an injured meniscus.

III. APPLICATION INSTRUCTIONS

- Make sure that the size of the device is properly selected.
- Before putting on the brace, unfasten and loosen all the straps, (1a).
- Place your leg in the middle of the open brace and adjust its shape to the shape of your leg, make sure that the patella cut-outs are in the right place, (2b).
- Fasten the straps. First, while holding the brace in place, fasten the small straps (1 and 2), which will prevent the brace from opening and moving. Then fasten straps 3/4/5/6, (3c).
- Tighten the straps enough to hold the brace in its correct position. Remember that the fastened straps should not cause any difficulties in blood flow. If necessary, adjust the fastening and tension of the straps. Figure 4 shows a correctly applied brace.

SIZE CHART

Size chart:			
Size	(A) Thigh circumference 15 cm above the center of the patella	(B) Calf circumference 15 cm below the center of the patella	
S	39–47 cm	30,5–35,5 cm	
M	47–53 cm	35,5–40,5 cm	
L	53–59,5 cm	40,5–47 cm	
XL	59,5–67 cm	47–50,5 cm	

IV. COMPOSITION:

- Poliuretan 62%,
- Nylon 32%,
- Poliester 6%.

V. TERMS OF USE

- The device may only be used as recommended by a physician,
- The method of putting on the brace should be demonstrated by qualified medical personnel, in the presence of the person taking care of the user (if applicable),
- The device and its components should be used only for their intended purpose,
- If pain, swelling, sensory changes or any unusual reactions occur during the use of the device, stop using it immediately and consult your physician,
- Use of the device at night is only permissible if recommended by the attending physician,
- Using the brace for a long period of time may result in limitation of active and passive movements in the joints and muscle atrophy, therefore it is recommended to systematically perform exercises selected by qualified medical personnel as recommended by the attending physician,
- The device should be kept clean following the recommendations given in these instructions for use and on the label placed on the device. Improper cleaning and maintenance of the device may lead to skin lesions,
- The device cannot be used if any part is missing. Two of the same or different devices cannot be used at the same time.

In case of any serious incident that has occurred in relation to the device, please report to the manufacturer and the competent authority of the Member State in which the user or patient is established.

VI. INDICATIONS

The device can be used in the event of:

- full postoperative immobilization,
- securing the knee joint in order to gradually restore mobility,
- patellar tendon injury,
- meniscus injury,
- patella fracture/tear.

VII. CONTRAINDICATIONS

- Deep vein thrombosis of the leg and thigh,
- Soft tissue defects and wound healing problems,
- Existing problems with lymph drainage,
- A too tight brace may lead to increased pain, while a too loose brace will not provide the intended effect,
- In the event of skin changes in the areas where the brace adheres, its use is limited and possible only under medical supervision,
- It is necessary for the attending physician to assess the use of the brace in people with memory disorders of various origins, people with mental disorders and children,
- If tingling, numbness or pain occurs while using the brace, the brace should be removed and the physician should be contacted immediately.

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It is recommended to wash the device regularly. Before washing, remove the metal parts of the device and fasten all hook and loop fasteners of the brace.

- Hand wash in water with a neutral, anti-allergic detergent, at 30 degrees,
- Rinse in lukewarm water, do not use fabric softeners,
- Do not wring,
- Do not bleach,
- Do not iron,
- Do not dry clean,
- Do not tumble dry,
- Dry flat in a shaded and airy place - do not expose the device to direct sunlight. The device is not heat resistant and loses its properties at high temperatures above 120°C.

IX. WARRANTY

The manufacturer assures the customer that the item sold is brand new and of high quality.
The device is covered by a warranty (detailed warranty conditions can be found on the Manufacturer's website).

X. SYMBOL DESCRIPTION

	Manufacturer		Keep dry
	CE marking		Fragile, handle with care
	Catalogue number		Do not cut
	Date of manufacture		Do not use if package is damaged and consult instructions for use
	Consult instructions for use		Recyclable packaging
	Caution		Keep away from Sunlight
	Batch code		Hand wash
	Medical device		Do not bleach
	Unique Device Identifier		Do not iron
	Quantity		Do not dry clean
	Model number		Dry vertically
	Date of issue		Do not tumble dry
	This way up		

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