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INSTRUCTIONS FOR USE
Hip orthosis with adjustable flexion and abduction angle „T-Max Hip” O-1370

I. INTRODUCTION

Dear Customers,
Thank you for choosing a product of REHA FUND Sp. z o.o., a leading supplier of medical supplies and rehabilitation products on the Polish market.
The instructions for use contain tips and suggestions, the knowledge of which is essential for the correct use of device.
The photos used in these instructions are for illustration purposes only!

NOTE! Read the instructions for use before using the product. The manufacturer shall not be held liable for improper use of the product, failure to observe the rules of safe use or use other than for its intended purpose.

II. SPECIFICATIONS

The hip orthosis with flexion and abduction angle adjustment „T-Max Hip” O-1370 consists of:

- rigid outer shell for optimum support;
- ankle pad, allowing limited flexion and extension movement, as well as abduction and adduction at the hip joint;
- adjustable straps for additional adjustment;
- highly permeable fabric and skin-friendly material to ensure comfort during use.

Medical device intended for use on the right and left leg.

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III. INTENDED USE AND APPLICATION

The product is intended for people with history of injuries such as fractures or sprains, during rehabilitation after surgical treatment of the hip joint, and as preventive care for osteoarthritis of the hip joint.

V. CONTRAINDICATIONS

- The product is not recommended for use in cases of risk of thrombosis, advanced varicose veins, peripheral arterial disease or significant deformities of the joints or surrounding tissues, when there are sensory abnormalities in the areas of contact between the orthosis and the skin, when several orthoses are used at the same time, when the orthosis is used in heavy physical work, when the elderly have fragile skin, diabetes or vascular diseases.
- **Sensory disturbances, e.g. tingling** – may occur when the orthosis has put direct pressure on the body.
- **Muscle atrophy** – may occur when movement of the legs is reduced, leading to muscle weakness and reduced muscle mass.
- **Circulatory system complications** – pain, cooling of the limb below the orthosis, cyanosis, varicose vein formation and superficial and deep vein phlebitis.

VI. ADVERSE REACTIONS

- **Mechanical abrasions and wounds in the area of contact between the orthosis and the skin** – may occur when the orthosis is not properly fitted, when the orthosis is worn for long periods of time in hot conditions, when there is a significant deformity of the joints or surrounding tissues, when there are sensory abnormalities in the areas of contact between the orthosis and the skin, when several orthoses are used at the same time, when the orthosis is used in heavy physical work, when the elderly have fragile skin, diabetes or vascular diseases.
- **Sensory disturbances, e.g. tingling** – may occur when the orthosis has put direct pressure on the body.
- **Muscle atrophy** – may occur when movement of the legs is reduced, leading to muscle weakness and reduced muscle mass.
- **Circulatory system complications** – pain, cooling of the limb below the orthosis, cyanosis, varicose vein formation and superficial and deep vein phlebitis.

VII. WARNINGS AND PRECAUTIONS

- The orthosis is intended for use by a single patient only. It must not be reapplied to another person.
- Consult your healthcare practitioner before and during use – especially if you have pain, swelling, skin changes or other disturbing symptoms.
- The device should only be used as directed by a doctor or medical personnel.
- Any modification of the device, including changing the angle of lower limb abduction is prohibited – all adjustments are made by the doctor.
- Do not attempt to repair the product or modify it in any way. Repairs and disassembly may only be carried out by authorized personnel.
- The product should be regularly inspected for wear and tear, damage and the condition of the surfaces in contact with the skin. In the event of damage, discontinue use immediately.
- Do not use the product on dry or damaged skin, open wounds or if you are allergic to any of the materials used.
- Do not use if you have sensory impairment, memory impairment, mental illness or other limitations – without the supervision of a responsible person.
- Remove the device during prolonged rest (e.g. sleep) or before a radiological examination.
- Do not use the product in a manner incompatible with its intended use.
- Store in a dry, clean, well-ventilated room, away from sources of heat, moisture and corrosive gases.
- Do not use accessories or add-ons not recommended by the manufacturer or a doctor.
- Contactant use with orthopedic devices is only possible after medical consultation.

VIII. GENERAL RULES FOR SAFE USE

Before use, consult a healthcare professional on the purpose of the application, the selection of the device and how to use it.

Each 'serious incident' related to the product shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient resides.

IX. METHOD OF USE

Make sure the product is not damaged before use. Once removed from the packaging, open the plastic buckles and Velcro to fit the orthosis onto the lower limb.

9.1 Instructions for fitting the orthosis

Fig. 1 Wrap the hip belt around the pelvis. The top edge of the belt should rest above the iliac crest.
Fig. 2 You may need help from other people from the top edge in Fig. 2. Wrap the thigh strap around the thigh. The lower edge should rest above the kneecap.
Fig. 3 Place the metal frame next to the pelvis and thigh. Locate the axis of the hinge over the greater trochanter (the most protruding bony point on the upper lateral side of the thigh). The top edge of the upper section should not be above the waist belt. The plastic plates of the lower frame should be within the thigh strap.




Fig. 4 Put the buckles on the belt and fit the straps tightly.
Fig. 5 Secure the thigh straps snugly around the thigh belt—tight enough to hold firmly, but not so tight as to cause discomfort.

9.2 Length adjustment instructions:

The rule applies to the upper and lower parts

Use the Allen key to adjust the flexion angle. Six motion lock positions limit hip flexion from 0° to 100° (0°, 20°, 40°, 60°, 80°, 100°).

Use an Allen key to adjust the abduction/adduction angle.

9.4 Instructions for adjusting the angle of abduction:

Use an Allen key to adjust the abduction/adduction angle.

X. CLEANING AND MAINTENANCE

The orthosis should be cleaned with mild detergents (water and soap). Do not chlorinate, do not chemically clean, do not wring out. Dry away from artificial heat sources. Dry undressed. Do not iron.

XI. STORAGE AND TRANSPORT

Store between 10°C and 25°C, 40-60% RH in a well ventilated place. Keep away from sharps to prevent damage.

Transport between 10°C and 25°C, 40-60% RH is possible by land, sea or air.

XII. END OF USE

If any abnormalities are found when inspecting the product in accordance with the manufacturer's recommendations see: Cleaning and Maintenance – step using the product immediately.

Disposal

At the end of use, the product must be disposed of in accordance with local environmental regulations.

XIII. INTENDED PERIOD OF USE

The intended lifetime of the product is 3 years
This time may be longer or shorter, depending on the frequency of use and maintenance.

XIV. WARRANTY

The manufacturer assures the customer that the product sold is brand new and of high quality. The product is covered by a warranty (see the manufacturer's website for detailed warranty conditions). The product should be used as intended and as recommended in the instructions for use.

Manufacturer	LOT	Batch code
		Quantity
		Unique product identification code
Catalogue No.		Temperature Limitation
Date of manufacture		Relative humidity limit
Consult instructions for use		Protect from humidity
		Hand wash
Medical device		Do not iron
Wash with soap		Do not bleach
		Do not dry-clean
		Dry horizontally
		Do not tumble dry

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